Pulmodyne®

BiTrac Select™

Interchangeable

Anti-Asphyxia

Bronchoscopy Elbow

22mm Female

EU*

Intersurgical UAB
Tel: +370 387 66 611
e-mail: info@intersurgical.lt
Amionijų g. 60, LT-18170, Pabrada, Lithuania

*European Union

Single Patient Use
Made in USA

Pulmodyne Inc.
2055 Executive Drive
Indianapolis, IN 46241
USA

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Intended Use

The Interchangeable Anti-Asphyxia Bronchoscopy Elbow is an accessory to the BiTrac Select Full Face Mask, which allows CPAP and Bilevel therapy to be delivered to the patient during a bronchoscopy procedure. Use of this elbow with these masks is an accessory for use with CPAP and Bilevel systems used in the treatment of adult Obstructive Sleep Apnea (OSA) and/or ventilator support. A minimum pressure of >3.0cm of H₂O at the mask is required.

Note

The Interchangeable Anti-Asphyxia Bronchoscopy Elbow is a temporary use elbow for use only during a bronchoscopy procedure.

Contraindications

1. This device should not be used on uncooperative or unresponsive patients.
2. This device may not be suitable for use on patients unable to remove the mask or are taking medication which may cause vomiting.

Warnings

1. This elbow should be used only during a Bronchoscopy procedure while in the presence of trained medical personnel.
2. This elbow does not have an integral leak port; therefore a separate exhalation device must be part of the circuit. Utilize the equipment manufacturer's recommended exhalation port in the designed circuit.
3. This device is designed to be used only in conjunction with a BiTrac Select Full Face Interchangeable Elbow compatible mask.

Before Use

Verify entrainment valve function on Interchangeable Anti-asphyxia Bronchoscopy Elbow.

- With the unit turned OFF; the entrainment valve flap should be open, allowing user to breathe room air
- Turn CPAP or Bilevel unit ON
- The flap valve should close, allowing the user to receive airflow from the unit
- If the valve does not function in this manner; replace the elbow
- Check valve before each use for malfunction or blockage caused by patient secretions

Use of Interchangeable Elbow Feature:

1. Remove the existing elbow from the mask by holding the adapter ring and twisting the elbow counter clockwise to disengage the bayonet lock (Refer to Figure B.)
2. Gently pull the elbow from the mask. (Refer to Figure C.)
3. Gently press the desired Interchangeable Compatible Elbow into the mask adapter ring and turn clockwise until fully seated. Do not use excessive force.

Bayonet Lock References:

Figure A.

Figure B.

Figure C.

Operating Specifications: 5°C to 40°C at humidity range of 15% to 95%

Storage Specifications: -20°C to 60°C at relative humidity up to 95% non-condensing

Caution: Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.

1 of 15

2691E-06-15-2022

Pulmodyne®

BiTrac Select™

Interchangeable

Anti-Asphyxia

Bronchoscopy Elbow

22mm Female

EU*

Intersurgical UAB
Tel: +370 387 66 611
e-mail: info@intersurgical.lt
Amionijų g. 60, LT-18170, Pabrada, Lithuania

*European Union

Single Patient Use
Made in USA

Pulmodyne Inc.
2055 Executive Drive
Indianapolis, IN 46241 USA

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