# Pulmodyne®

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**da** Brugsanvisning

**no** Bruksanvisning

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Made in Malaysia

Made in Malaysia

AccuPAP is a Trademark of Pulmodyne®

\*European Union

Caution: Federal Law (USA) restricts this device to sale by  
or on the order of a licensed physician.

2891E-06-21-2022

# AccuPAP™ Directions for Use

## DESCRIPTION:

The AccuPAP Positive Airway Pressure therapy system is a non-sterile Single Patient Use, respiratory therapy device. The AccuPAP Positive Airway Pressure therapy system has the ability to provide supplemental oxygen when supplied by compressed oxygen.

The AccuPAP system consists of:

1. Positive Airway Pressure Device
  - Ambient air entrainment orifice
  - Supply gas inlet connection
  - Pressure monitoring port with cap
2. Patient interface - mouthpiece
3. 7 foot supply tubing for connecting the AccuPAP to the source gas.

## INDICATIONS FOR USE:

The AccuPAP is a Positive Airway Pressure therapy system that is indicated for the prevention or reversal of atelectasis. The AccuPAP Positive Airway Pressure therapy system is for patients who are capable of following directions. It may also be used to optimize delivery of aerosolized medication when utilized with a small volume nebulizer.

## CONTRAINDICATIONS:

Although no absolute contraindications to the use of PAP Therapy have been reported, the following should be carefully evaluated before a decision is made to initiate therapy:

- Inability to tolerate increased work of breathing
- Intracranial pressure (ICP) >20mm Hg
- Recent facial, oral or skull surgery or trauma
- Esophageal surgery
- Untreated pneumothorax
- Known or suspected tympanic membrane rupture or other middle ear pathology
- Hemodynamic instability
- Acute sinusitis
- Epistaxis
- Active hemoptysis
- Nausea

## WARNINGS:

1. Use of this device at excessive pressures may have adverse effects. Expert clinical judgment should be exercised in the selection of settings for each individual patient.
2. **DO NOT OCCLUDE AMBIENT AIR ENTRAINMENT ORIFICE.**
3. **DO NOT EXCEED 15 LPM GAS INPUT.**

## Adverse reactions may include:

Increased work of breathing that may lead to hypoventilation and hypercarbia; increased cranial pressure; cardiovascular compromise; myocardial ischemia; decreased venous return; air swallowing with increased risk of vomiting and aspiration; claustrophobia; skin break down and discomfort if used in conjunction with a mask; and pulmonary barotrauma.

## PRECAUTIONS:

1. DO NOT AUTOCLAVE.
2. DO NOT IMMERSE COMPLETELY IN ANY LIQUID.
3. DO NOT DISASSEMBLE EXCEPT FOR MOUTHPIECE OR TUBING.

## INSTRUCTIONS FOR USE:

### 1. Set-up:

- a. Connect one end of tubing (supplied with the device) to a 0-15 LPM flow meter that is connected to a 50-60 psi regulated medical air or oxygen (if prescribed) supply source.
- b. Connect the other end of the tubing to the device's supply gas inlet connection.
- c. If pressure monitoring is not being done, seal port with attached pressure monitoring port cap.
- d. Attach the mouthpiece to the rounded patient interface end of the device.

### 2. Initial Flow Meter Setting:

- a. Instruct the patient to relax while performing diaphragmatic breathing.
- b. Set initial flow rate to 8 LPM on the air or oxygen supply source flow meter. Adjust flow meter to desired level as needed.
- c. Set the desired airway pressure to the desired expiratory airway pressure.
- d. Place mouthpiece in patient's mouth and instruct patient to breathe easily against the pressure from the device.
- e. Slow inhalations and exhalations by the patient will maintain adequate airway pressure during the breathing cycle.
- f. Maintain breathing through the AccuPAP Positive Airway Pressure Therapy System to accomplish clinical goals.
- g. Several repetitions may be needed to ensure that individual patient needs are being met.

### 3. AccuPAP Procedure:

- a. Ensure the flow meter is set as identified in previous steps.
- b. If tolerated, have patient sit with elbows resting comfortably on table.
- c. Place mouthpiece lightly into mouth.
  - Be sure patient maintains a tight seal on the mouthpiece during therapy.
  - Recommend the use of a nose clip, if necessary.
- d. Instruct the patient to breathe diaphragmatically.
- e. Perform breathing through the device as recommended. AccuPAP Positive Airway Pressure Therapy System is intended to be used on an intermittent basis (not for prolonged periods of time).
- f. Repeat steps b to e as prescribed.

## FOR USE WITH DISPOSABLE PRESSURE MANOMETER (if provided)

### Instructions:

1. Attach the manometer with its flexible connector onto the device's pressure monitoring port.
2. Occlude the patient port on the attached device and pressurize the manometer to the maximum scale reading on the manometer.
3. Release the pressure and check that the float returns smoothly back to "0" cm H2O mark.
4. Attach the system to the patient and monitor the patient to assure proper ventilation.
5. For a point of reference, set the "O" ring to be approximately at the desired pressure reading.
6. When repeated use is required for the same patient, check the function and accuracy of the manometer prior to each use.
7. Discard after use.

### WARNINGS:

- For Single Patient Use. Do not clean or sterilize as this may affect accuracy or function of the manometer.
- Check function and accuracy prior to each use, including when repeated use is required for the same patient.
- A minimal amount of leakage of airway gas is normal. The effects of this leakage must be evaluated for each patient.

### CAUTIONS:

- The device must be used only by personnel thoroughly trained in the techniques of pulmonary resuscitation or airway management.
- During use, the patient's condition must be monitored. Should the performance of this device appear to be erratic during use, the device must be removed and replaced as required.

### Manometer Accuracy:

- ±1 cm H2O from 0-10 cm H2O
- ±2 cm H2O from 10-40 cm H2O
- ±3 cm H2O from above 40 cm H2O

## FOR USE WITH NEBULIZER

Suggested Instructions For Use With Small Volume Nebulizer:

### CAUTION:

It is important that device instructions for the AccuPAP Positive Airway Pressure Therapy System, as well as the intended nebulizer, are understood prior to use. The selected nebulizer should be used according to the manufacturer's instruction for use in a positive pressure system.

**DO NOT OCCLUDE.** Inspect device on a routine basis to ensure proper use and function. If damaged, do not use. Verify all connections are secure.

### Set-up:

Follow set-up instructions for each device.

### DISPOSAL:

Dispose of AccuPAP Positive Airway Pressure Therapy System in a safe manner according to Federal/State/Local regulations and guidelines for disposal of contaminated medical waste.

### Ordering Information:

Each system includes AccuPAP, one 7ft. connecting tube, mouthpiece, one pressure port cap and if desired, optional disposable manometer: one of the following:

AccuPAP system without disposable manometer (10 per Case) .....313-6000  
AccuPAP system with disposable manometer (10 per Case) .....313-6001

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