# Pulmonary<sup>®</sup>

## O2-MAX™ Trio System

Pulmonary's O2-MAX System delivers continuous positive airway pressure (CPAP) throughout the breathing cycle. It provides CPAP at preset levels throughout inspiration and exhalation, independent of the patient's flow rate. O2-MAX System is intended for use on spontaneously breathing patients, there is no assembly required. **SINGLE PATIENT USE ONLY.**

### INDICATIONS FOR USE:

To provide CPAP to spontaneously breathing adult (>30kg) patients in the hospital and pre-hospital (EMS) environment.

### CONTRAINDICATIONS:

May be contraindicated for patients with any of these conditions:

- Facial lacerations
- Laryngeal trauma
- Recent tracheal or esophageal anastomosis
- Gastrointestinal bleeding or ileus
- Recent gastric surgery
- Basilar skull fracture
- Patients at high risk of vomiting
- Emphysematous Bulla - when an area of the lung may be brittle and present a risk of bursting
- Hypovolaemia - low blood volume

**Operating Specifications:** 5°C to 40°C at humidity range of 15% to 95%

**Storage Specifications:** -20°C to 60°C at relative humidity up to 95% non-condensing

### HOW IT WORKS:

O2-MAX Trio Generator is a variable flow venturi device that uses an oxygen supply in conjunction with entrained air to generate an output flow. O2-MAX Generator uses a 4.5 bar oxygen supply, and can generate flows up to 140 lpm and fractional inspired oxygen (FiO<sub>2</sub>) at levels of approximately 30%, 60%, 90%+. The O2-CPAP valves, which are snapped onto the anti-asphyxia housing end of the circuit, are used to maintain positive pressure flow rates from 60 to 140 lpm.

### OPERATING THE SYSTEM:

- Connect directly to a 4.5 bar gas source
- Prior to use, be sure the device is free of obstructions and verify proper valve function
- The Trio generator device is supplied in the 30% FiO<sub>2</sub> position. To adjust for different FiO<sub>2</sub> levels, turn the knob until desired level is visible in the window. The printed percentage position must be fully visible in order to achieve the desired level. If not fully visible, approximately 30% will be delivered.
- Monitor with any in-line pressure gauge or manometer with a range of 0-30cm H<sub>2</sub>O per local established protocol. If the pressure drops significantly during patient inspiration then the flow is too low, increase the flow rate to the patient.
- For use only by thoroughly trained personnel
- Do not use O2-CPAP valve if it becomes occluded. An occluded O2-CPAP valve may obstruct patient's exhalation and result in potential injury. If this occurs, discard entire system or remove the occluded O2-CPAP valve from the system and replace with another O2-CPAP valve.
- The circuit has an anti-asphyxia valve which prevents exhalation into the tubing in the event of no fresh gas flow. If the fresh gas flow fails, the anti-asphyxia valve in the circuit will minimize the risk of asphyxia when there is no oxygen pressure. When there is not oxygen flow from the generator, the mask should not be worn.

### WARNING:

- Do not connect any gas supply other than oxygen to the O2-MAX System
- When not in use, disconnect O2-MAX oxygen line from 4.5 bar source.
- Watch the O2-CPAP valve to ensure that it remains open during inspiration
- Monitor with any in-line pressure gauge or manometer with a range of 0-30cm H<sub>2</sub>O per local established protocol. If the pressure drops significantly during patient inspiration then the flow is too low, increase the flow rate to the patient.
- For use only by thoroughly trained personnel
- Do not use O2-CPAP valve if it becomes occluded. An occluded O2-CPAP valve may obstruct patient's exhalation and result in potential injury. If this occurs, discard entire system or remove the occluded O2-CPAP valve from the system and replace with another O2-CPAP valve.
- The circuit has an anti-asphyxia valve which prevents exhalation into the tubing in the event of no fresh gas flow. If the fresh gas flow fails, the anti-asphyxia valve in the circuit will minimize the risk of asphyxia when there is no oxygen pressure. When there is not oxygen flow from the generator, the mask should not be worn.

### MONITORING THE PATIENT:

During operation, be sure to check the following on a regular basis:


- Ensure that there are no leaks at the patient connection
- Ensure that there is flow from the O2-CPAP valve during inspiration (which means that the generator is supplying adequate flow to meet patient demand). Optionally, monitor the in-line pressure gauge or manometer during inspiration. If the pressure drops then the flow is inadequate
- Monitor the patient's arterial blood oxygen saturation (SaO<sub>2</sub>)
- Monitor the patient's SpO<sub>2</sub> continuously using pulse oximetry
- Monitor the patient for signs of dehydration and discomfort in the upper airways

### CAUTIONS:

Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.

Single Patient Use Only. Do not sterilize or immerse the O2-MAX System or any of its components in any solution. Dispose of O2-MAX System according to local established protocols when finished with use.

<p><b>System includes:</b></p> <p>Trio Generator with Filter 72" Anti-Asphyxia Circuit O2-CPAP™ Valve Mask with Head Strap</p>		<p><b>Pulmonary Inc.</b> 2055 Executive Drive Indianapolis, IN 46241 USA</p> <p><a href="http://www.pulmonary.com">www.pulmonary.com</a></p> <p>EC REP QNET BV Kantstraat 19 NL 5016 NP Haaren The Netherlands</p>
--	---	--

\*European Union

Single Patient Use

Made in the USA

Pulmonary and ...*bringing change to life* are Registered Trademarks of Pulmonary, Inc. O2-MAX, Trio and O2-CPAP are trademarks of Pulmonary, Inc.

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