

System includes:

1. Fixed Flow Generator with Filter
2. 72" Anti-Asphyxia Circuit
3. O2-CPAP™ Valve
4. Mask with Head Strap
5. O2-MAX (Guard)

1. Place the smooth side of the O2-MAX Guard Filter onto the anti-asphyxia housing of the circuit.
2. Snap into place (Ref. A)
3. Ensure the PEEP Valve is fully seated onto the filter before beginning therapy (Ref. B)
4. Check for airflow and function as part of circuit checkout procedure, prior to use.

WARNING: Replace filter at least every 24 hours, or earlier if increased resistance is noted.

CAUTIONS:

Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.

Single Patient Use Only. Do not sterilize or immerse the O2-MAX System or any of its components in any solution. Dispose of O2-MAX System according to local established protocols when finished with use.

Pulmodyne Inc.
2055 Executive Drive
Indianapolis, IN 46241 USA

KNET BV
Kantstraat 19
NL 5076 NP Haaren
The Netherlands

Intersurgical UAB
Tel: +370 897 66611
e-mail: info@intersurgical.lt
Amšionų g. 60, LT-18170, Pakruoje, Lithuania

Medical Device

Latex Free


Made in the USA

Single Patient Use

Pulmodyne and ...bringing change to life are Registered Trademarks of Pulmodyne, Inc. O2-MAX and O2-CPAP are trademarks of Pulmodyne, Inc.

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Operating Specifications: 5 °C to 40 °C at humidity range of 15% to 95%

Storage Specifications: -20 °C to 60 °C at relative humidity up to 95% non-condensing

O2-MAX Generator is a fixed flow venturi device that uses an oxygen supply in conjunction with air to generate an output flow. O2-MAX Generator uses a 50psi (~4 bar) oxygen supply, and can generate flows up to 140 lpm and fractional inspired oxygen (FiO₂) approximately 30%. The O2-CPAP valves, which are snapped onto the anti-asphyxia housing end of the circuit, are used to maintain positive pressure flow rates from 60 to 140 lpm.

Pulmodyne's O2-MAX Fixed Generator can be converted to a variable flow unit when used in conjunction with the O2-MAX Stand Alone Variable Control Unit (sold separately).

OPERATING THE SYSTEM:

- Connect directly to a 50psi (~4 bar) gas source.
- Prior to use, be sure the device is free of obstructions and verify proper valve function.
- Place mask over patients face. Utilize the head strap to secure the mask firmly in place.

WARNING:

- Do not connect any gas supply other than oxygen to the O2-MAX System.
- When not in use, disconnect O2-MAX oxygen line from 50psi (~4 bar) source.
- Watch the O2-CPAP valve to ensure that it remains open during inspiration.
- Monitor with any in-line pressure gauge or manometer with a range of 0-30cm H₂O per local established protocol. If the pressure drops significantly during patient inspiration then the flow is too low, increase the flow rate to the patient.
- For use only by thoroughly trained personnel.
- Do not use O2-CPAP valve if it becomes occluded. An occluded O2-CPAP valve may obstruct patient's exhalation and result in potential injury. If this occurs, discard entire system or remove the occluded O2-CPAP valve from the system and replace with another O2-CPAP valve.
- The circuit has an anti-asphyxia valve which prevents exhalation into the tubing in the event of no fresh gas flow. If the fresh gas flow fails, the anti-asphyxia valve in the circuit will minimize the risk of asphyxia when there is no oxygen pressure. When there is not oxygen flow from the generator, the mask should not be worn.

MONITORING THE PATIENT:

During operation, be sure to check the following on a regular basis:

- Ensure that there are no leaks at the patient connection.
- Ensure that there is flow from the O2-CPAP valve during inspiration (which means that the generator is supplying adequate flow to meet patient demand). Optionally, monitor the in-line pressure gauge or manometer during inspiration. If the pressure drops then the flow is inadequate.
- Monitor the patient's arterial blood oxygen saturation (SaO₂).
- Monitor the patient for signs of dehydration and discomfort in the upper airways.

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