## Laryngo-Tracheal Mucosal Atomization Device

**INDICATIONS FOR USE:** for atomizing solutions across the naso, laryngeal, and oropharyngeal mucous membranes.

### **DIRECTIONS FOR USE:**

1. Fill syringe with desired volume of solution and eliminate remaining air.
2. Connect the DART-Reach device to the syringe.
3. Place the tip of the DART in the nostril, laryngeal, or oropharyngeal cavity.
4. Compress the syringe plunger to spray atomized solution into the nasal, laryngeal, or oropharyngeal cavity.
5. Repeat steps if needed.

**WARNING:** The use of products intended for single patient use on more than one patient may result in degraded performance or loss of functionality and may also result in exposure to viral, bacterial, fungal, or prionic pathogens. Validated cleaning and sterilization methods and instructions for reprocessing to original specifications are not available for these products. This product is not designed to be cleaned, disinfected or sterilized.

<b>SPECIFICATIONS:</b>	
AVERAGE PARTICLE SIZE (MMAD)	50µm
APPLICATOR DEAD SPACE	REACH600/700 = 0.22mL REACH720 = 0.13mL
TIP DIAMETER	0.19"
APPLICATOR LENGTH	REACH600/700 = 8" REACH720 = 4"
SYRINGE ACCURACY (REACH600)	± 5%

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