Intended use:

This Depth Filter is designed for use with ventilators, anesthesia machines and open flow systems where filtration of inspired and/or expired gases is desired. Product Specifications and Information are listed in the table below.

Set Up:

1. Place the Filter in the circuit in the desired location, at the machine end, connected to the Inspiratory outlet and/or expiratory inlet, or between the artificial airway and the proximal breathing circuit.
2. Attach gas sampling line to luer port, if present.
3. Ensure connections are secure.
4. Check for airflow and function as part of circuit checkout procedure prior to use.

Cautions:

1. Do not resterilize, soak, rinse or reuse.
2. Ensure all connections are secure at all times.
3. Replace unit immediately if there is any contamination, occlusion or any indication of malfunction.
4. The dead space of this product should be taken into consideration when determining tidal volume and patient ventilation requirements.
5. Dispose of properly.
6. Federal (USA) law restricts this device to sale by or on the order of a physician.

Warnings:

1. Replace Filter at least every 24 hours or earlier if increased resistance is noted.

Contraindications:

1. Filter should not be used in the proximal airway with patients producing fulminating frothy secretions within their airways, or patients with hemoptysis.
2. Do not use in conjunction with conventional humidifiers.
3. Do not add moisture to the Filter.
4. During delivery of inhaled medication the Filter should be removed or bypassed.

Part Number Sterile/Clean	Description	Resistance to Flow (cmH2O @ L/min)	Dead Space (mL)	Weight (g)	Bacterial Efficiency	Viral Efficiency	Tidal Volume (mL)
200201 / 200202	Aero-Pro™ Filter Small Angled	1.4 @ 30	21.1	15.6	99.99	99.95	150-1000
200601 / 200602	Aero-Pro™ Filter Small Straight	1.4 @ 30	19.9	14.4	99.99	99.95	150-1000
200401 / 200402	Aero-Pro™ Filter Compact Angled	2.0 @ 60	35.4	23.5	99.999	99.99	150-1000
200301 / 200302	Aero-Pro™ Filter Compact Straight	2.0 @ 60	31.5	23.7	99.999	99.99	150-1000
200101 / 200102	Aero-Pro S™ Filter	1.0 @ 60	46.9	29.5	99.998	99.97	N/A
313-435	Bacterial/Viral Depth Filter No Port	1.0 @ 60	30.0	22.5	99.9998	99.9998	>125

Note: No natural rubber or synthetic latex was used in the manufacturing of this product.

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