# Pulmodyne®

## BiTrac SE™

### Nasal Mask with Standard Elbow 22mm Female

#### Intended Use

The Pulmodyne BiTrac Nasal mask and accessories (mask, forehead pad and head strap) are used with and attached to a positive pressure air source, CPAP or Bilevel equipment, and has the Indications for Use consistent with the Indications for Use of the equipment to which it is attached.

#### Warnings

- 1) This mask does not have an integral leak port; therefore a separate exhalation device must be part of the circuit. Utilize the equipment manufacturer's recommended exhalation port in the designated circuit.
- 2) At low CPAP Pressures the flow through the exhalation port may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- 3) If skin irritation develops, this may be due to excessive pressure or skin reaction or a combination of the two. Loosen the head strap and/or try another mask or forehead pad size.
- 4) This CPAP mask should only be used with CPAP systems recommended by your physician or respiratory therapist. A mask should not be used unless the CPAP system is turned ON and operating properly. The vent hole or holes associated with the mask should never be blocked.

*Explanation of the Warning: CPAP systems are intended to be used with the special masks with connectors which have vent holes to allow continuous flow of air out of the mask. When the CPAP machine is turned on and functioning properly, new air from the CPAP machine flushes the exhaled air out through the attached mask exhalation port. However, when the CPAP machine is not operating, enough fresh air will not be provided through the mask, and exhaled air may be rebreathed. Rebreathing of exhaled air for longer than several minutes can in some circumstances lead to suffocation. This warning applies to most models of CPAP systems.*


#### Pulmodyne

2055 Executive Drive  
Indianapolis, IN 46241 USA

www.Pulmodyne.com

EC REP QNET BV  
Kantstraat 19  
NL-5076 NP Haaren  
The Netherlands

Intersurgical UAB  
Tel: + 370 387 66 611  
e-mail: info@intersurgical.lt  
Armoniu g. 60, LT-18170, Pabradė, Lithuania

\*European Union

Latex Free

#### Single Patient Use

Made in USA

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.

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