BiTrac™ Circuit Non-Vented

Bilevel Application

Caution: Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.

Latex Free

Single Patient Use

Pulmodyne Inc.
2055 Executive Drive
Indianapolis, IN 46241
USA

T: +1.317.246.5505
F: +1.317.246.5501

www.pulmodyne.com

Pulmodyne is a registered trademark of Pulmodyne, Inc. BiTrac and Bilevel Compatible are trademarks of Pulmodyne, Inc.

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Made in
USA

72" Smooth ID Tube with attached filter, swivel connector

Intended Use

The BiTrac Circuit is used with and attached to a positive pressure air source, CPAP or Bilevel equipment, and has the Indications for Use consistent with the Indications for Use of the equipment to which it is attached.

Warnings

1. This circuit does not have an integral leakport; therefore a separate exhalation device must be part of the circuit/patient interface.
2. Before instituting mechanical ventilation, confirm the proper connection of the line circuit and verify ventilator operation. Patient, ventilator and circuit must be monitored on a regular basis according to local medical protocol.
3. Do not pull or stretch the circuit tubing. This could result in circuit leakage.
4. This circuit is not recommended for use with a heated humidifier since it does not allow for use of water traps.

This circuit is designed for Single Patient Use. It is not intended to be cleaned or sterilized. The length of time a circuit is used should be determined by local medical protocol.

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