Pulmodyne® O2-RESQ™ Circuit

Pulmodyne's O2-RESQ Circuit delivers continuous positive airway pressure (CPAP) throughout the breathing cycle. It provides CPAP at preset levels throughout inspiration and exhalation, independent of the patient's flow rate. O2-RESQ Circuit is intended for use on spontaneously breathing patients.

INDICATIONS FOR USE:

To provide CPAP to spontaneously breathing adult (>30kg) patients in the hospital and pre-hospital (EMS) environment.

CONTRAINDICATIONS:

May be contraindicated for patients with any of these conditions:

- Facial lacerations
- Laryngeal trauma
- Recent tracheal or esophageal anastomosis
- Gastrointestinal bleeding or ileus
- Recent gastric surgery
- Basilar skull fracture
- Patients at high risk of vomiting
- Emphysematous Bulla - when an area of the lung may be brittle and present a risk of bursting
- Hypovolaemia - low blood volume

Operating Specifications: 5°C to 40°C at humidity range of 15% to 95%

Storage Specifications: -20°C to 60°C at relative humidity up to 95% non-condensing

HOW IT WORKS:

The O2-RESQ circuit is intended for use with the MyO CPAP-system. MyO's operation is based on the Venturi-principle, using an oxygen supply in conjunction with entrained air to generate an output flow. The MyO uses a 46 to 85 psi oxygen supply, can generate flows up to 140 lpm and allows for fractional inspired oxygen (FiO2) to be set between 30% and nearly 100%. The O2-CPAP valves, which are snapped onto the anti-asphyxia housing end of the circuit, are used to maintain preset positive pressure at flow rates from 60 to 140 lpm.

OPERATING THE SYSTEM:

- Connect the MyO-system to the oxygen supply, then listen for leaks
- Connect the filter to the Venturi-port of the MyO
- Connect the corrugated tubing to the patient port of the MyO and the pressure monitoring tubing to the MyO pressure port
- Prior to use, check to be sure the device is free of obstructions and verify proper valve function
- Set the O2-CPAP valve according to the patient's needs
- Place mask over patient's face. Utilize the head strap to secure the mask firmly in place

WARNING:

- Do not connect any gas supply other than oxygen to the O2-RESQ circuit.
- Watch the O2-CPAP valve to ensure that it remains open during inspiration.
- Check the MyO-monitor for changes in pressure. If the pressure drops significantly during patient inspiration, then the flow is too low - increase the flow rate to the patient.
- For use only by thoroughly trained personnel.
- Do not use O2-CPAP valve if it becomes occluded. An occluded O2-CPAP valve may obstruct patient's exhalation and result in potential injury. If this occurs, discard entire system or remove the occluded O2-CPAP valve from the system and replace it with another O2-CPAP valve.
- The circuit has an anti-asphyxia valve which prevents exhalation into the tubing in the event of no fresh gas flow. If the fresh gas flow fails, the anti-asphyxia valve in the circuit will minimize the risk of asphyxia when there is no oxygen pressure. When there is no oxygen flow from the generator, the mask should not be worn.

MONITORING THE PATIENT:

During operation, be sure to check the following on a regular basis:

- Ensure that there are no leaks at the patient connection
- Ensure that there is flow from the preset O2-CPAP valve during inspiration (which means that the MyO-generator is supplying adequate flow to meet patient demand). Optionally, check the MyO-monitor for changes in pressure. If the pressure drops then the flow is inadequate.
- Monitor the patient's arterial blood oxygen saturation (SaO2)
- Monitor the patient for signs of dehydration and discomfort in the upper airways.
- Using the MyO, monitor the patient's delivered FiO2


CAUTIONS:

Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.

Single Patient Use Only. Do not sterilize or immerse the O2-RESQ Circuit or any of its components in any solution. Dispose of O2-RESQ Circuit according to local established protocols when finished with use.

Circuit includes:

1. 72" Corrugated Anti-Asphyxia Circuit with accompanying pressure monitoring tubing
2. O2-CPAP™ Valve
3. Mask with Head Strap
4. Filter

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Pulmodyne
2055 Executive Drive
Indianapolis, IN 46241 USA

www.Pulmodyne.com

*European Union

Intersurgical UAB
Tel: + 370 387 66611
e-mail: info@intersurgical.lt
Amioniu g. 60, LT-18170, Pabradė, Lithuania

Single Patient Use

Pulmodyne is a Registered Trademark of Pulmodyne, Inc.
O2-RESQ and O2-CPAP are trademarks of Pulmodyne, Inc.

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