## BiTrac SE MaxShield™ Shield Mask with Standard Elbow 22mm Female

Single Patient Use

Made in USA

\*European Union

EU\*

Intersurgical UAB  
Tel: + 370 387 66 611  
e-mail: info@intersurgical.lt  
Armonijų g. 60, LT-18170, Pabradė, Lithuania

Latex Free

**Pulmodyne Inc.**

2055 Executive Drive  
Indianapolis, IN 46241 USA

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Caution: Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.

Pulmodyne and BiTrac SE Max Shield are trademarks of Pulmodyne, Inc.

### Intended Use

The Pulmodyne BiTrac SE Shield mask is an accessory for use with CPAP or Bilevel Systems for ventilatory support.

### Contraindications

- 1) This device should not be used on uncooperative or unresponsive patients.
- 2) This device may not be suitable for use on patients unable to remove the mask or are taking medication which may cause vomiting, have nocturnal vomiting, hiatal hernia or impaired cardiac sphincter function.
- 3) This mask may not be suitable for use on patients with the following conditions: glaucoma, recent eye surgery or dry eyes.

### Warnings

- 1) Use only with a ventilator equipped with adequate alarms and safety systems for ventilator failure. This mask does not include an anti-asphyxia valve that would allow a patient to breathe if the ventilator failed.
- 2) This mask does not have an integral leak port; therefore a separate exhalation device must be part of the circuit. Utilize the equipment manufacturer's recommended exhalation port in the designated circuit.
- 3) This device is not suitable for providing life support ventilation.
- 4) Any unusual chest discomfort, shortness of breath, stomach distension, belching, severe headache, drying eyes, eye pain, eye infections or blurred vision when receiving positive airway pressure or upon awakening should be reported to the clinician.
- 5) To minimize risks of vomiting during sleep, avoid eating or drinking for three (3) hours before use.
- 6) Make sure the mask and headgear are the correct size for the patient's face. If excessive discomfort is experienced, adjust the mask and head gear. If the discomfort continues, consider an alternative mask.

**Operating Specifications:**  
5°C to 40°C at humidity range  
of 15% to 95%

**Storage Specifications:**  
-20°C to 60°C at relative  
humidity up to 95% non-  
condensing

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