


AccuPAP™

Positive Airway Pressure Therapy System

The AccuPAP system provides PEEP and supplemental oxygen in a handheld device. The integrated PEEP valve provides a true mechanical peep to support lung recruitment. It may also be used with a T-piece to optimize delivery of aerosolized medication when utilized with a small volume nebulizer. The AccuPAP is a great therapy to treat patients with atelectasis, or for preventative care.



Positive Airway Pressure Device

PEEP valve (5/10/15/20cm H₂O)

Flow set at 8LPM



DPM™
Disposable
Pressure
Manometer

AccuPAP - Positive Airway Pressure Therapy System

Part No.	Description	Case/Qty.
313-6000	AccuPAP with Mouthpiece and 7' Oxygen Tubing	10
313-6001	AccuPAP with Mouthpiece, 7' Oxygen Tubing and DPM™	25