

# DART™

## Intranasal Mucosal Atomization Device

Deliver medication quickly and painlessly with the DART. Atomized medication delivery results in therapeutic onset delivery similar to intravenous delivery. There's no special training required for a system wide adoption. The delivery method is painless for the patient, helping them remain calm in a high stress situation.



*ALS and BLS approved*

*Measure and deliver with accuracy*

*Patient comfort*



**DART100**



**DART110**



**DART300**

### DART-Intranasal Mucosal Atomization Device

Part No.	Description	Case/Qty.
DART100	DART Nasal Device with 3mL Syringe	25
DART110	DART Nasal Device with 1mL Syringe	25
DART140	DART Nasal Device with 3mL Syringe and Vial Adapter	25
DART300	DART Nasal Device without Syringe	25



**DART140**  
\*Vial Adapter  
(Blunt Plastic Needle)

4255E-03-10-2022