

NASO-FLO®

The Naso-Flo® nasopharyngeal airway device allows for direct inline oxygen delivery to the upper pharynx. The humidification vents positioned towards the distal tip facilitate heat and moisture transfer. It can also be supplied with an optional respiratory indicator and integrated capnography.










Integrated Capnography

Adjustable size

Direct inline oxygen delivery

Integrated Capnography line



4611B-03-10-2022

Pulmodyne®

www.pulmodyne.com | 317.900.4093

Twitter Facebook LinkedIn Instagram
@pulmodyne

Naso-flo® Nasopharyngeal Airway Device							
Part No.	O2-Port	15mm Connector	Size	Respiratory Indicator	Orange Connector	Adjustable	Case/Qty.
313-NF414	Yes	Yes	ID 4.0 / 18 FR	Yes	Yes	Yes	10
313-NF415	Yes	Yes	ID 5.0 / 22 FR	Yes	Yes	Yes	10
313-NF416	Yes	Yes	ID 6.0 / 25 FR	Yes	Yes	Yes	10
313-NF417	Yes	Yes	ID 7.0 / 29 FR	Yes	Yes	Yes	10
313-NF404	Yes	Yes	ID 4.0 / 18 FR	Yes	Yes	-	10
313-NF405	Yes	Yes	ID 5.0 / 22 FR	Yes	Yes	-	10
313-NF406	Yes	Yes	ID 6.0 / 25 FR	Yes	Yes	-	10
313-NF407	Yes	Yes	ID 7.0 / 29 FR	Yes	Yes	-	10
313-NF408	Yes	Yes	ID 8.0 / 32 FR	Yes	Yes	-	10
313-NF409	Yes	Yes	ID 9.0 / 36 FR	Yes	Yes	-	10
313-NF204	Yes	Yes	ID 4.0 / 18 FR	Yes	-	-	10
313-NF205	Yes	Yes	ID 5.0 / 22 FR	Yes	-	-	10
313-NF206	Yes	Yes	ID 6.0 / 25 FR	Yes	-	-	10
313-NF207	Yes	Yes	ID 7.0 / 29 FR	Yes	-	-	10
313-NF208	Yes	Yes	ID 8.0 / 32 FR	Yes	-	-	10
313-NF209	Yes	Yes	ID 9.0 / 36 FR	Yes	-	-	10
313-NF104	Yes	Yes	ID 4.0 / 18 FR	-	-	-	10
313-NF105	Yes	Yes	ID 5.0 / 22 FR	-	-	-	10
313-NF106	Yes	Yes	ID 6.0 / 25 FR	-	-	-	10
313-NF107	Yes	Yes	ID 7.0 / 29 FR	-	-	-	10
313-NF108	Yes	Yes	ID 8.0 / 32 FR	-	-	-	10
313-NF109	Yes	Yes	ID 9.0 / 36 FR	-	-	-	10