

EZ-Spray™


The EZ-Spray was designed for rapid anesthetization of patients. It is a single patient use, closed system, preventing the possibility of cross-contamination. Adjust the oxygen source to operate at 7-10 LPM according to desired flow rate. The 180° directional tip provides the ability to effortlessly direct and deliver medication to the upper pharynx. It is a quick and thorough anesthetization process.



Disposable & Single Patient Use

Power Spray Atomizer

Cross-Contamination Prevention



EZ-Spray - Powered Atomization Device

| Part No. | Description | Case/Qty. |
|----------|---------------------------------------|-----------|
| 313-100 | EZ-Spray (60cc) with 7' Oxygen Tubing | 10 |
| 313-101 | EZ-Spray (15cc) with 7' Oxygen Tubing | 10 |
| 313-102 | EZ-Spray SMR (15cc) | 10 |

3636B-08-02-2022