

Quik-Connects™

Disposable Oxygen Connections

Connect to an oxygen source instantly with the Quik-Connects. These single-use disposable DISS gas fittings are designed to be used with the O2-MAX to prevent cross-contamination. They are available in 9 different connections, including the United States and International. These disposable, single-use products can quickly transport a patient from bedside to the ambulance and between hospital rooms.



Disposable & Single-Use


JCAHO Standard

Cross-contamination prevention



Quik-Connect Disposable DISS Gas Fittings

Part No.	Description	Case/Qty.
313-111111	OHMEDA-MALE-DISS	30
313-121111	CHEMTRON-MALE-DISS	30
313-131111	PURITAN BENNETT-MALE-DISS	30
313-141111	SCHRADER-MALE-DISS	30
313-351111	AMERICAN GAS-MALE-DISS	30
313-361111	GERMAN-MALE-DISS	30
313-371111	AFNOR-MALE-DISS	30
313-381111	BRITISH-MALE-DISS	30
313-391111	UNI-MALE-DISS	30



Compatible with the O2-MAX™

8291E-03-10-2022